




## 3 STAGE BALLOON DILATOR OTW

Intended for endoscopic dilation of airway tree. Supplied sterile in a peel open pouch. Designed for single use only.



### Key features:

- One balloon provides three different pre - defined diameters – 6-7-8mm, 8-9-10mm, 10-11-12mm, 12-13.5-15mm, 15-16.5-18mm
- Non compliant material offers optimum radial force and uniform inflated balloon diameter.
- Radiopaque bands to identify balloon location fluoroscopically
- Atraumatic tip prevents damage to the endoscope
- Accepts 0.035" guidewire

### Specifications

Size

**Sheath Diameter** 6.8 Fr

**Length** 1800 and 2400 mm

**Balloon Length** 55 mm