




HF ELECTRODE - BALL ELECTRODE 3mm and 5 mm

3 mm and 5 mm Ball Electrode is used to coagulate and smoothen the freshly resected surface of the prostate gland during TURP procedures. Product supplied Sterile and single use only.



Key features:

- High quality PTFE insulation.
- Designed to be used with monopolar current.
- Each electrode individually tested for High Frequency leak test.
- Compatible with standard resectoscopes.

Classification:

Compatible with

- STORZ 1 Stem, 30 Degrees Instrument
- STORZ 2 Stem, 30 Degrees Instrument
- WOLF 1 Stem, 15 Degrees Instrument

- OLYMPUS 1 Stem, 12 Degrees Instrument
- OLYMPUS 1 Stem, 30 Degrees Instrument

Other properties:

Available as Disposable loop as well as Multiuse loop

Specifications

Size

Diameter 24 Fr