




## **BILIARY PLASTIC STENTS**

Ace Biliary Stents are used to drain obstructed Biliary ducts. This device is supplied sterile in a peel-open pouch. Intended for single use only.



### **Key features:**

- Stents made of polyethylene and / or PTFE
- Smooth surface for less trauma to the tissue.
- Highly Radiopaque for better visualization.
- Rounded Taper ends.
- Flap and Pigtail design for better retention.
- High Flow Rate

### **Classification:**

- Standard
- Amsterdam

- Straight
- Tennabaum
- Double Pigtail
- Single Pigtail
- Pancreatic

## Specifications

### Size

**Diameter** 5, 7, 8.5, 9,10, 11.5 & 12 Fr

**Length** 3 - 20 cm