




## HF ELECTRODE - COLLINS KNIFE

Collins Knife is used to make straight incisions in a very hard prostate tissue during TURP procedures. Product supplied Sterile and single use only.



### Key features:

- High quality extruded PTFE insulation.
- Each electrode individually tested for High Frequency leak test.
- Designed to be used with monopolar current.
- Compatible with standard resectoscopes

### Classification:

Compatible with

- STORZ 1 Stem, 30 Degrees Instrument
- STORZ 2 Stem, 30 Degrees Instrument
- WOLF 1 Stem, 25 Degrees Instrument
- OLYMPUS 1 Stem, 12 Degrees Instrument
- OLYMPUS 1 Stem, 30 Degrees Instrument

## Other properties:

Available as Disposable loop as well as Multiuse loop

## Specifications

Size

**Diameter** 24 Fr