




## HF ELECTRODE - CUTTING LOOP

Cutting Loops are used to remove prostate tissue in the form of small chips during TURP procedures. Product supplied Sterile and single use only.



### Key features:

- High quality PTFE insulation.
- Designed to be used with monopolar current.
- Each electrode individually tested for High Frequency leak test.
- Compatible with standard resectoscopes.

### Classification:

Compatible with

- STORZ 1 Stem, 30 Degrees Instrument
- STORZ 2 Stem, 30 Degrees Instrument
- WOLF 1 Stem, 25 Degrees Instrument
- OLYMPUS 1 Stem, 12 Degrees Instrument
- OLYMPUS 1 Stem, 30 Degrees Instrument

## Other properties:

Available as Disposable loop as well as Multiuse loop.

## Specifications

Size

**Diameter** 24 Fr