




GUIDE WIRES - PTFE

AcePTFE Guidewires are used for percutaneous vessel entry to facilitate subsequent introduction(s) of an intravascular device. This device is Supplied sterile in a peel-open pouch. Intended for single use only..



Key features:

- Fixed core, Flexible floppy tip - for atraumatic advancement through urinary tract
- Stainless Steel with PTFE jacket - easy tracking of stents and urological instruments
- Excellent torqueability for precise tip control
- Low of coefficient of friction for smooth advancement
- Better pushability for increased placement control
- Extremely smooth.
- Safe vessel access.

Classification:

Available in straight and J Tip variants

Specifications

Size

Diameter 0.018", 0.025", 0.028", 0.032", 0.035" & 0.038"

Length 150 cm

Other

Balloon Max Fluid Injection