



LEEP/LLETZ LOOP ELECTRODES - BALL

The ACE LEEP/LLETZ electrodes are intended for use in the excision, coagulation or electro-conization of abnormal tissue growth or lesions in the cervix called cervical intraepithelial neoplasia (CIN) and HPV related condylomata acuminata. The Ball Electrode is used for Fulguration and Desiccation. Product supplied Sterile and single use only.



Key features:

- Standard 3/32" (2.3 mm) shaft
- Stable at High temperatures.
- For Single Use Only.

Specifications

Size

Diameter 3, 4 and 5 mm

Shaft Diameter 2.3 mm

Length 13 cm