



LAPAROSCOPIC MONOPOLAR J HOOK - DISPOSABLE

The Ace disposable Laparoscopic J Hook is intended for use in minimally invasive laparoscopic surgical procedures for coagulating, cauterizing and other manipulation of tissues and vessels during laparoscopic procedures.



Key features:

- High Quality Insulation.Ensures insulation down to the blades.
- Dielectric strength 20 KV/mm.
- Ergonomic Molded Handle.
- **India's FIRST and ONLY Manufacturer of Single Use Instruments.**

Classification:

Available Tip Length 3mm and 5 mm

Specifications

Size

Diameter 5 mm

Length 330 mm

Tip Length 3 ~ 5 mm