



ACE MULTIPLE BAND LIGATOR SET / KIT

The ACE Multiple Band Ligator is a variceal ligation system used for the endoscopic ligation of esophageal varices. Intended for single use only.



Key features:

- **Ace Multi Band Ligator** offers a consistent high quality and is a user friendly system
- **Ace Multi Band Ligator** provides an exact performance substitute to the best in the market
- Least compression of the endoscope as minimum release force needed for bands
- Latex Free Bands
- 4,5 or 6 high quality bands

Classification:

The product is available in various sizes to fit endoscopes ranging from 8.5 mm to 11 mm diameter. Supplied

un – Sterile.

AMBLS-6 Set with 6 bands for endoscope diameter 8.5 ~ 11.0 mm

AMBLS-5 Set with 5 bands for endoscope diameter 8.5 ~ 11.0 mm

AMBLS-4 Set with 4 bands for endoscope diameter 8.5 ~ 11.0 mm

For Sterile Sets: add ‘S’ to catalogue number.

Specifications

Size

Diameter 8.5 - 11 mm