




## BILIARY METALLIC STENTS

**Used for Stenosis or strictures in the Biliary tract. Supplied sterile in a peel open pouch with deployer.**



### **Key features:**

- Available as Fully Covered, Both Ends Bare, With removal string
- Available TTS and as Percutaneous transhepatic
- High corrosion resistance and shape memory
- Slow expansion below body temperature
- Unique deploying mechanism is reconstrainable over 98% of deployment to aid in accurate repositioning of the stent
- Optimum Radial Force
- Atraumatic stent endings to reduce tissue trauma
- Multiple radiopaque markers to facilitate accurate stent placement

## Specifications

### Size

**Diameter** 6,8,10 mm

**Length** 40, 50, 60, 70, 80,100,120 mm