



## ENDOPYELOTOMY STENTS

Used for temporary internal drainage from the ureteropelvic junction (UPJ) to the bladder following incision of a stricture.



### Key features:

- Segment with no sideports prevents post-op narrowing of the ureteral lumen
- Prevents growth of ureteral wall to the stent
- Highly radiopaque for better visualization
- Large diameter at incised Utero Pelvic Junction provides optimal healing & drainage
- Pigtail shape for better retention.

### Specifications

#### Size

**Diameter** 4/8 , 5/10, 6/10, 6/12, 6/14, 7/14, 7/16 Fr

**Length** 16 - 24 - 26 - 28 cm