




## ENUCLEATION BUTTON ELECTRODE - BIPOLAR

Ace Enucleation Button Electrode is used for bipolar enucleation, coagulation and vaporization for treatment of benign prostate hyperplasia (BPH). Compatible with Wolf and Storz 1 Stem Instrument. Product supplied Sterile and single use only.



### Key features:

- Fast bipolar enucleation, coagulation and vaporization.
- Minimal penetration of thermal energy.
- Improved visibility and enhanced durability during procedures.

### Classification:

Compatible with:

Storz 1 Stem Instrument

Wolf 1 Stem Instrument

## Specifications

### Size

**Diameter** 22, 24 Fr