




## EUS - FNB NEEDLE

An Endoscopic Ultra Sound (EUS) - Fine Needle Biopsy (FNB) is used to sample targeted submucosal Gastrointestinal lesions through the accessory channel of an ultrasound endoscope. This device is supplied sterile and is intended for single use only. To be used by or on order of a physician only.



### Key features:

- Enhanced Echogenic needle tip for optimal ultrasonic visualization.
- Unique needle design for higher and consistent tissue yield.
- Ergonomically contoured handle for better control and stability.
- Actuating lock for safer depth control of needle.
- Ultra smooth sheath with adjustable length attuned for a full range of endoscopes.

### Classification:

## EUS Needle Ordering Information

Needle Type	Needle Size Ø	19G	22G	25G
ACE SUPER CORE		AEUS-SFNB-19	AEUS-SFNB-22	AEUS-SFNB-25
ACE MASTER CORE	FNB	AEUS-MFNB-19	AEUS-MFNB-22	AEUS-MFNB-25
ACE ULTRA CORE		AEUS-UFNB-19	AEUS-UFNB-22	AEUS-UFNB-25
ACE DUAL CORE		AEUS-DFNB-19	AEUS-DFNB-22	AEUS-DFNB-25
ACE SUPER DUAL CORE	FNB	AEUS-SDFNB-19	AEUS-SDFNB-22	AEUS-SDFNB-25
ACE MASTER DUAL CORE		AEUS-MDFNB-19	AEUS-MDFNB-22	AEUS-MDFNB-25
ACE ULTRA DUAL CORE		AEUS-UDFNB-19	AEUS-UDFNB-22	AEUS-UDFNB-25

## Specifications

### Size

**Sheath Diameter** 5.2 Fr

**Needle Size** 19, 22 and 25G