



## GUIDE WIRES - STRIPED

Ace Striped Guidewires are used for percutaneous vessel entry to facilitate subsequent introduction(s) of an intravascular device. This device is Supplied sterile in a peel-open pouch. Intended for single use only..



### Key features:

- Straight, Flexible tip - for atraumatic advancement through urinary tract
- Nitinol core with striped jacket - easy tracking of stents and urological instruments
- Low coefficient of friction for smooth advancement
- Better pushability for increased placement control
- Extremely smooth.
- Safe vessel access.

### Classification:

with / without hydrophilic floppy tip

### Specifications

#### Size

**Diameter** 0.018", 0.025", 0.028", 0.032", 0.035" & 0.038"

**Length** 150 cm