



HEMORRHOIDAL BAND LIGATOR

Used to ligate internal hemorrhoids, facilitated by an anoscope. Supplied unsterile. For single use only.



Key features:

- Anoscope with in-built LED light
- Anoscope available in two sizes – Medium and Large
- Supplied with 4 latex free rubber bands
- Suction control at finger tip
- Automatic Band counter

Classification:

CATALOGUE NUMBERS	DEVICE DESCRIPTION
AHBL-04	ACE Hemorrhoidal Band Ligator 04 Bands Without Anoscope

AHBL-04AM	ACE Hemorrhoidal Band Ligator 04 Bands With Medium Anoscope
AHBL-04AL	ACE Hemorrhoidal Band Ligator 04 Bands With Large Anoscope