



## HYDROPHILIC GASTRO GUIDEWIRE

Used to facilitate the placement of instruments during diagnostic or treatment procedures. Supplied sterile in a peel open Tyvek pouch. For single use only



### Key features:

- Unmatched torqueability
- Excellent trackability.
- Extremely smooth
- Safe vessel access
- Nitinol non kink core

### Specifications

#### Size

**Diameter** 0.018, 0.021, 0.025, 0.028, 0.032, 0.035 Inch

**Length** 260 and 450 cm