



LOOP BASKETS

Used for endoscopic removal of foreign bodies, stones, food boluses and severed tissues such as polyps.. Supplied sterile for single use only.



Key features:

- Reinforced nylon net weave offers durability, strength and visibility throughout the procedure
- Designed for smooth and concise removal of foreign bodies even in torturous positions
- Highly effective for removing coins, watch batteries or other foreign bodies that are commonly swallowed by mistake
- Ergonomic handle that can be rotated 360°

Specifications

Size

Diameter 7 Fr

Basket Diameter 20,30,40 mm

Length 230 cm

