



NASAL BILIARY DRAINAGE CATHETER (NBDC)

Used for temporary drainage of the Biliary duct through the use of an indwelling catheter. Supplied sterile in a peel open pouch. Designed for single use only.



Key features:

- Pigtail shape for better retention.
- Each set includes nasal transfer tube and a connector to facilitate irrigation and drainage
- Available as Straight, Straight with Pigtail and Angled with Pigtail variants

Classification:

Specifications

Size

Diameter 5, 6, 7, 8, 8.5 & 10 Fr

Length 250 cm