




POLYPECTOMY SNARES

Used for Endoscopic Excision of polyps using a monopolar electrosurgical unit and retrieval of foreign bodies from gastrointestinal tract. Supplied sterile for single use only.



Key features:

- Available in Hexagonal, Oval and Crescent shapes
- Designed for smooth and concise polyp excision even in difficult positions
- Ergonomic handle that can be rotated 360°
- Available in a variety of configurations to accommodate procedural requirements

Specifications

Size

Diameter 6, 13, 15, 20, 25, 30, 35 mm

Length 230 cm