



# SPHINCTEROTOME / PAPILLOTOME

Used for Endoscopic Cannulation of Biliary Ductal System and Sphincterotomy.



## Key features:

- Triple Lumen with pre curved tapered Tip
- Cutting Wire Length: 20, 30 mm
- Tip design: 5.5 Fr
- Available with / without insulation
- Efficient, smooth and easy to use handle
- Recommended Guidewire 0.035 inch.

## Specifications

### Size

**Sheath Diameter** 6.8 Fr

**Length** 1800 and 2000 mm