



## STONE RETRIEVAL BASKET NITINOL

**Intended to retrieve Foreign Bodies and C.B.D. stones endoscopically. Designed for single use only.**



### **Key features:**

- Engineered to facilitate pushability
- Unique handle for ease of use
- Atraumatic tip for ease of procedure
- Compatible with 2.8 mm instrument channel
- Hexagonal, Trapezoidal and Balloon shaped

### **Specifications**

#### **Size**

**Diameter** 15, 20, 25 30 and 35 mm

**Length** 2300 mm