




## TRI PRONG GRASPER

Used to facilitate endoscopic manipulation/ removal of stones and retrieval of foreign objects. Designed for single use only.



### Key features:

- Designed for minimum trauma.
- Maximum grasping efficiency.

### Classification:

For URS and Cystoscope

- Rigid
- Flexible
- PCNL

## Other properties:

Rigid

Size : 1.9, 2.5, 3, 4, 5 and 6

Length : 40 and 60 cms

Flexible

Size : 1.9, 2.5, 3, 4, 5 and 6

Length : 40 and 60 cms

PCNL

Size : 7, 9 and 12 Fr

Length : 40 cms

## Specifications

Size

**Diameter** 1.9 - 12 Fr

**Length** 40 - 60 cm