



URETERAL ACCESS SHEATH

Used to establish a conduit during endoscopic urological procedures facilitating the passage of flexible scopes into the urinary tract. Designed for single use only.



Key features:

- Atraumatic tip for easy access
- Kink resistant coil for smooth navigation through a torturous anatomy
- Irradium-Platinum marker for highest radiopacity
- Tapered dilator for easy intramural access
- Can be supplied with hydrophilic coating

Specifications

Size

Diameter 9.5/11.5, 10.5/12.5, 12/14 Fr

Length 35, 45 and 55 cm