



## URETERAL ACCESS SHEATH

Used for ureteral dilation and to establish a conduit during endoscopic urological procedures facilitating the passage of flexible scopes into the urinary tract. Designed for single use only



### Key features:

- Atraumatic tip for easy access
- Kink resistant coil for smooth navigation through a torturous anatomy
- Iridium-Platinum marker for highest radiopacity
- Tapered dilator for easy intramural access
- Can be supplied with hydrophilic coating

### Specifications

#### Size

**Diameter** 9.5/11.5, 10.5/12.5, 12/14 Fr

**Length** 28, 35, 38, 45 and 55 cm

