


Ace Stone Retrieval Devices

Helical Baskets

Designed for the manipulation and retrieval of stones from the Bladder, Ureter and Renal Pelvis. This device is supplied sterile in a peel-open pouch. Intended for one-time use.



Cat. No.	Size	Length (x)
H – 3-4-x	3 Fr, 4 wire	70, 90, 110 cm
H – 4-4-x	4 Fr, 4 wire	70, 90, 110 cm
H – 4-5-x	4 Fr, 5 wire	70, 90, 110 cm
H – 5-4-x	5 Fr, 4 wire	70, 90, 110 cm
H – 5-5-x	5 Fr, 5 wire	70, 90, 110 cm


X is the length of the Stone Retrieval.
For Stainless Steel add SS to Cat. No.
For Nitinol Baskets add N to Cat. No.



Ace Stone Retrieval Devices

Flat Wire Baskets

Designed for the manipulation and retrieval of stones from the Bladder, Ureter and Renal Pelvis. This device is supplied sterile in a peel-open pouch. Intended for one-time use.



Cat. No.	Size	Length (x)
S - 3-4-x	3 Fr, 4 wire	70, 90, 110 cm
S - 4-4-x	4 Fr, 4 wire	70, 90, 110 cm
S - 4-5-x	4 Fr, 5 wire	70, 90, 110 cm
S - 5-4-x	5 Fr, 4 wire	70, 90, 110 cm
S - 5-5-x	5 Fr, 5 wire	70, 90, 110 cm


X is the length of the Stone Retrieval.
For Stainless Steel add SS to Cat. No.
For Nitinol Baskets add N to Cat. No.



Ace Stone Retrieval Devices

Tipless Baskets

Designed to make stone retrieval easy in very difficult areas even at the base of the kidney's calyces. It protects against mucosal damage and perforation. This device is supplied sterile in a peel-open pouch. Intended for one-time use.



Cat. No.	Size	Length (x)
T – 3-4-x	3 Fr, 4 wire	70, 90, 110 cm
T – 4-4-x	4 Fr, 4 wire	70, 90, 110 cm
T – 4-5-x	4 Fr, 5 wire	70, 90, 110 cm
T – 5-4-x	5 Fr, 4 wire	70, 90, 110 cm
T – 5-5-x	5 Fr, 5 wire	70, 90, 110 cm

X is the length of the Stone Retrieval,

For Stainless Steel add SS to Cat. No. For Nitinol Baskets add N to Cat. No.

For Helical shape add H to Cat. No. For Flat Wire shape add S to Cat. No.



Ace Stone Retrieval Devices

Tri Prong Grasper

Designed to facilitate endoscopic manipulation and /or removal of stones as well as for retrieval of foreign objects. Tri prong forceps are designed to be less traumatic while providing maximum grasping efficiency. This device is supplied sterile in a peel-open pouch. Intended for one-time use.



Cat. No.	Size	Length (x)
GF – 3-x	3 Fr Tri Prong Grasping Forceps	70, 90, 110 cm
GF – 4-x	4 Fr Tri Prong Grasping Forceps	70, 90, 110 cm
GF – 5-x	5 Fr Tri Prong Grasping Forceps	70, 90, 110 cm

X is the length of the Grasper

