

ACE URETERAL DJ STENTS


Ureteral Stents (Closed Distal Tip)



Used for temporary internal drainage from the Uretropelvic junction to the bladder. Closed distal tip for back loaded placement. This device is supplied sterile in a peel-open pouch. Intended for one-time use.

- Highly radiopaque for better visualization and has an extremely low encrustation tendency.
- Soft polyurethane material ensures patient comfort.
- Wide lumen ensures a superior drainage.
- Stents are tapered at both ends, that helps in easy penetration to the exact location.
- Can be kept indwelling for 120 days. However periodic evaluation is strongly advised.
- **Carbothane Long Life Stents** are also available in white and blue colour.

Caution : Periodic Evaluation is advised.



| Sr No | Size in French | Length in cm | Catalogue Nos (DJ size/length) |
|-------|----------------|--------------|--------------------------------|
| 01. | 3 Fr | 08 – 30 cm | DJ 3/ (length) |
| 02. | 4 Fr | 08 – 30 cm | DJ 4/ (length) |
| 03. | 4.5 Fr | 08 – 30 cm | DJ 4.5/ (length) |
| 04. | 4.7 Fr | 08 – 30 cm | DJ 4.7/ (length) |
| 05. | 5 Fr | 08 – 30 cm | DJ 5/ (length) |
| 06. | 5.5 Fr | 08 – 30 cm | DJ 5.5/ (length) |
| 07. | 6 Fr | 08 – 30 cm | DJ 6/ (length) |
| 08. | 6.5 Fr | 08 – 30 cm | DJ 6.5/ (length) |
| 09. | 7 Fr | 08 – 30 cm | DJ 7/ (length) |
| 10. | 8 Fr | 08 – 30 cm | DJ 8/ (length) |

* Please specify length * Each Stent is supplied with a Pusher

For Ureteral Stents (Closed End) with pusher, with Guidewire : (DJ size/length) G

For Ureteral Stents (Closed End) with pusher, with Suture : (DJ size/length) S

For Ureteral Stents (Closed End) with pusher, with Suture & Guidewire : (DJ size/length) GS

For Ureteral Stents (Closed End) with pusher, with Clamp & Guidewire & Suture: (DJ size/length) GSC

For Ureteral Stents (Closed End) Multiloop One End : (DJO size/ML)

For Ureteral Stents (Closed End) Multiloop Both End : (DJO size/ML B)


Ureteral Stents (Open Distal Tip)



Used for temporary internal drainage from the Uretropelvic junction to the bladder. Open taper tip facilitates easier over the wire placement. This device is supplied sterile in a peel-open pouch. Intended for one-time use.

- Highly radiopaque for better visualization and has an extremely low encrustation tendency.
- Soft polyurethane material ensures patient comfort.
- Wide lumen ensures a superior drainage.
- Stents are tapered at both ends, that helps in easy penetration to the exact location.
- Can be kept indwelling for 120 days. However periodic evaluation is strongly advised.
- **Carbothane Long Life Stents** are also available in white and blue colour..

Caution : Periodic Evaluation is advised.



| Sr No | Size in French | Length in cm | Catalogue Nos (DJO size/length) |
|-------|----------------|--------------|---------------------------------|
| 01. | 3 Fr | 08 – 30 cm | DJO 3/ (length) |
| 02. | 4 Fr | 08 – 30 cm | DJO 4/ (length) |
| 03. | 4.5 Fr | 08 – 30 cm | DJO 4.5/ (length) |
| 04. | 4.7 Fr | 08 – 30 cm | DJO 4.7/ (length) |
| 05. | 5 Fr | 08 – 30 cm | DJO 5/ (length) |
| 06. | 5.5 Fr | 08 – 30 cm | DJO 5.5/ (length) |
| 07. | 6 Fr | 08 – 30 cm | DJO 6/ (length) |
| 08. | 6.5 Fr | 08 – 30 cm | DJO 6.5/ (length) |
| 09. | 7 Fr | 08 – 30 cm | DJO 7/ (length) |
| 10. | 8 Fr | 08 – 30 cm | DJO 8/ (length) |

* Please specify length * Each Stent is supplied with a Pusher

For Ureteral Stents (Closed End) with pusher, with Guidewire : (DJO size/length) G

For Ureteral Stents (Closed End) with pusher, with Suture : (DJO size/length) S

For Ureteral Stents (Closed End) with pusher, with Suture & Guidewire : (DJO size/length) GS

For Ureteral Stents (Closed End) with pusher, with Clamp & Guidewire & Suture: (DJO size/length) GSC

For Ureteral Stents (Closed End) Multiloop One End : (DJO size/ML)

For Ureteral Stents (Closed End) Multiloop Both End : (DJO size/ML B)

Mono J Urinary Diversion Stents

The Mono-J Urinary Diversion Stent is used for interoperative placement to stent the ureter during percutaneous, endoscopic or operative procedures. This device is supplied sterile in a peel-open pouch. Intended for one-time use.



| Sr No | Size in French | Length in cm | Guidewire | Open Tip | Closed Tip |
|-------|----------------|--------------|-----------|------------|------------|
| 01. | 4.7 | 70 | 0.032'' | MJO 4.7/70 | MJ 4.7/70 |
| 02. | 5 | 70 | 0.035'' | MJO 5/70 | MJ 5/70 |
| 03. | 6 | 70 | 0.038'' | MJO 6/70 | MJ 6/70 |
| 04. | 7 | 70 | 0.038'' | MJO 7/70 | MJ 7/70 |

* Each Stent is supplied with a Pusher


For Mono J Stents with Guidewire add suffix G to Cat Nos

Endopyelotomy Stents

Used for temporary internal drainage from the Uretropelvic junction to the bladder following incision of a stricture. The segment with no sideports prevents postoperative narrowing of the ureteral lumen while preventing ingrowth of the ureteral wall to the stent. This device is supplied sterile in a peel-open pouch. Intended for one-time use.

Caution : Periodic Evaluation is advised.






| Sr No | Stent Segment Size in French | Pigtail Stent Segment Size in French | Length in cm | Catalogue Nos (EPI size/length) |
|-------|------------------------------|--------------------------------------|--------------|---------------------------------|
| 01. | 10 | 05 | 22 - 30 | EPI 10/5 (Length) |
| 02. | 12 | 06 | 22 - 30 | EPI 12/6 (Length) |
| 03. | 14 | 07 | 22 - 30 | EPI 14/7 (Length) |
| 04. | 16 | 08 | 22 - 30 | EPI 16/8 (Length) |

* Please specify length * Each Stent is supplied with a Pusher

Tumor Stents

Used for temporary as well as for palliative treatments of internal urinary drainage when tumors are formed in the urinary tract. This device is supplied sterile in a peel-open pouch. Intended for one-time use.

Caution : Periodic Evaluation is advised.



| Sr No | Size in French | Length in cm | Guidewire | Catalogue Nos (TS size/length) |
|-------|----------------|--------------|-----------|--------------------------------|
| 01. | 05 | 22 – 30 cm | 0.035” | TS 5/ (Length) |
| 02. | 06 | 22 – 30 cm | 0.038” | TS 6/ (Length) |
| 03. | 07 | 22 – 30 cm | 0.038” | TS 7/ (Length) |

* Please specify length

For Tumor Stents with Guidewire add suffix G to Cat Nos

