


BIOPSY FORCEPS

Used to collect tissue endoscopically for histologic examination and/or for foreign body retrieval. Supplied Sterile in a peel open pouch.

The Ace Disposable Biopsy Forceps is intended to take tissue sample for Biopsy from the Gastrointestinal tract through a gastrointestinal endoscope.



Key features:

- India's first and only manufacturer of high-quality Disposable Biopsy Forceps.
- Easy opening and closing of the biopsy forceps even within the angled endoscope
- Minimized friction for trouble – free passage through the working channel of endoscope
- Smooth sheath with perfect strength and flexibility for advancement and control during scope passage through a torturous anatomy
- Clean and Precise tissue cut – no tearing, no crushing of tissues.
- Quality, reliability and peace of mind
- Available as Serrated (Alligator), Oval Cup, Standard types with and without Spike
- Available in multiple lengths to accommodate various procedural requirements

Specifications

Size

Diameter 1.8, 2.3 mm

Length 70, 100, 150, 170, 210, 230, 250 cm