




Inflation Device

Used for controlled balloon dilation, to inflate and deflate the balloon and to measure the pressure within the balloon during the procedure. Designed for single use only.



Key features:

- 30 ATM Pressure Gauge
- 20 ml clear transparent syringe
- Easy to operate lock – unlock mechanism
- Stable, Accurate pressures with a high degree of repeatability
- Efficient Intra-cam locking mechanism allows maintaining pressures and guarantees fast automatic release

Classification:

- With/without Three way stopcock
- With/without Hemostasis – Y connector kit.

