


ACE LAPROSCOPIC GRASPERs - DISPOSABLE

The Ace Indian Laparoscopic Monopolar Serrated Grasper is intended for use in Laparoscopic procedures involving grasping and holding.



Key features:

- High Carbon / Corrosion Resistance Stainless Steel Blades
- Dual- Edge Blades with a smooth cutting action
- Wide opening to take in larger chunks of tissue for cutting
- Does not loose its sharpness over extensive use
- Allows Dissection and transection of fine tissue
- High Quality Insulation Tubing Ensuring insulation down to the blades (Dielectric Strength - 20 Kv/mm)
- Ergonomic Molded Handle that is comfortable for any size of hand.
- Molded with a Banana Monopolar Electrosurgery Connector
- Designed for a Comfortable 360 degree Rotation by your Fingertip.
- *India's FIRST and ONLY Manufacturer of Single Use Instruments.*

Specifications

Size

Diameter 5 mm

Working Length 330 ~ 430 mm

Jaw Length 18 mm