


Logo ace medical devices

LEEP/LLETZ LOOP ELECTRODES - ROUND

The ACE LEEP/LLETZ electrodes are intended for use in the excision, coagulation or electro-conization of abnormal tissue growth or lesions in the cervix called cervical intraepithelial neoplasia (CIN) and HPV related condylomata acuminata. Product supplied Sterile and single use only.



4BDAB243 BBB4 C89B DEB26933987D7D4B

Key features:

- Standard 3/32" (2.3 mm) shaft
- Tungsten wire ensures an accurate and smooth excision..
- Maintains better shape and wire integrity.
- Stable at High temperatures.
- For Single Use Only.

Specifications

Size

Length 13 cm

Width 1.0 ~ 2.5 cm