

# LEEP/LLETZ LOOP ELECTRODES - ROUND

The ACE LEEP/LLETZ electrodes are intended for use in the excision, coagulation or electro-conization of abnormal tissue growth or lesions in the cervix called cervical intraepithelial neoplasia (CIN) and HPV related condylomata acuminata. Product supplied Sterile and single use only.



## Key features:

- Standard 3/32" (2.3 mm) shaft
- Tungsten wire ensures an accurate and smooth excision..
- Maintains better shape and wire integrity.
- Stable at High temperatures.
- For Single Use Only.

## Specifications

### Size

**Length** 13 cm

**Width** 1.0 ~ 2.5

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email [info@acemedicaldevices.com](mailto:info@acemedicaldevices.com) [acemedicaldevices.com](http://acemedicaldevices.com)

cm