


## **SURGICAL ADHESIVE**

Ace Surgical Adhesive is a sterile, ready-to-use cyanoacrylate adhesive to be used on the skin. It is a blue liquid compound that polymerizes rapidly in the presence of moisture and is used for the closure of uncomplicated skin lacerations. It is a surgical adhesive that provides a safe option for wound closure.



### **Key features:**

- Instantaneous solidification in contact with blood
- Forms a strong and flexible bond with tissue
- Negligible histotoxicity & non-carcinogenic
- Bacteriostatic and Biodegradable
- 0.5 ml closes wound upto 10 cms
- 1.0 ml closes wound upto 20 cms

### **Specifications**

#### **Size**

**Diameter**

**Length** 0.5 - 1 ml