


GUIDE WIRES - HYDROPHILIC

Ace Hydrophilic Guide wires and are used for percutaneous vessel entry to facilitate subsequent introduction(s) of an intravascular device. This device is Supplied sterile in a peel-open pouch. Intended for single use only.



Key features:

- Unmatched torqueability
- Excellent trackability.
- Extremely smooth
- Safe vessel access
- Nitinol non kink core

Classification:

Straight and Angled Tip

Specifications

Size

Diameter 0.018", 0.025", 0.028", 0.032", 0.035" & 0.038"

Length 150 cm