


## AMPLATZ RENAL DILATOR SET

Used for progressive dilatation for the nephrostomy tract prior to percutaneous kidney stone removal and to introduce the Amplatz sheath which is intended for percutaneous dilation of tissue surrounding a previously placed dilator.



### Key features:

- Radiopaque dilators & sheaths
- Short tapered tips
- Smooth surface reduces patient's trauma
- Introducer Catheter acts as a guide for dilators 14F-30F
- Introducer Catheter Sheath allows safety wire placement
- Mini Amplatz Set 6 – 22 Fr is available

### Classification:

Can also supply along with Guidewire and / or Puncture needle

## Specifications

### Size

**Diameter** 6 - 30 Fr

**Length** 20 - 30 cm