


BALL ELECTRODE - BIPOLAR

Ball Electrode is used to coagulate and smoothen the freshly resected surface of the prostate gland during TURP procedures, using bipolar current. Product supplied Sterile and single use only.



Key features:

- Faster firing, longer lasting
- High quality insulation.
- Specially designed for use with bipolar current.
- .Compatible with all standard resectoscopes.

Classification:

Compatible with

- STORZ 1 Stem Instrument
- STORZ 2 Stem Instrument
- Olympus 1 Stem Instrument
- Wolf 1 Stem Instrument

Other properties:

Available with / without bipolar cable.

Specifications

Size

Diameter 24 Fr