

# **BALL ELECTRODE - BIPOLAR**

Ball Electrode is used to coagulate and smoothen the freshly resected surface of the prostate gland during TURP procedures, using bipolar current. Product supplied Sterile and single use only.



## **Key features:**

- Faster firing, longer lasting
- High quality insulation.
- Specially designed for use with bipolar current.
- .Compatible with all standard resectoscopes.

## **Classification:**

Compatible with

- STORZ 1 Stem Instrument
- STORZ 2 Stem Instrument
- Olympus 1 Stem Instrument
- Wolf 1 Stem Instrument

## Other properties:

Available with / without bipolar cable.

## Specifications

### Size

**Diameter** 24  
Fr