


## URETERAL BALLOON DILATOR

Used for transluminal dilation of ureteral strictures or ureteral dilation prior to ureteroscopy or stone manipulation. Designed for single use only.



### Key features:

- Balloon dilation is a non surgical technique that expands the blocked urinary tract
- Radiopaque markers are designed on the catheter to be seen under fluoroscopy
- Tri fold balloon with low profile can be used for maximum no of procedures
- Short tapered distal tip allows close approach to ureteral stone.
- Inflation pressure upto 16 ATM
- Provides effective and uniform dilation

### Specifications

#### Size

**Shaft Diameter** 3.5 and 4.8 Fr

**Dilator Length** 40, 60 cm

**Balloon Length** 20, 40, 60 mm

**Width** 5, 6 mm