

URETERAL BALLOON DILATOR

Used for transluminal dilation of ureteral strictures or ureteral dilation prior to ureteroscopy or stone manipulation. Designed for single use only.



Key features:

- Balloon dilation is a non surgical technique that expands the blocked urinary tract
- Radiopaque markers are designed on the catheter to be seen under fluoroscopy
- Tri fold balloon with low profile can be used for maximum no of procedures
- Short tapered distal tip allows close approach to ureteral stone.
- Inflation pressure upto 16 ATM
- Provides effective and uniform dilation

Specifications

Size

Shaft Diameter 3.5 and 4.8 Fr

Dilator Length 40, 60 cm

Balloon Length 20, 40, 60

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@ccomedicaldevices.com ccomedicaldevices.com

mm

Width 5, 6 mm