


BIOPSY GUN

Used for taking biopsies in CT, Ultrasound, surgery, urology, nephrology and other internal medicine applications.



Key features:

- Disposable Semi automatic Biopsy needle
- Light weight and compact: spring loaded biopsy needle easier to maneuver
- One hand operation: easily operated with one hand
- No risk of losing specimen: once the cannula is fired the Specimen is locked into the window of the stylet assuring an intact core sample for diagnosis.
- Sharp beveled point allows quick easy penetration with less trauma to organs and surrounding tissues.

Classification:

14, 16 G x 10, 15 Cms - For Breast Biopsy
16, 18 G x 15, 20 Cms - For Liver & Kidney

18, 20 G x 15, 20 Cms - For Lung
18 G x 20, 25, 30 Cms - For Prostate

Specifications

Size

Diameter 14 - 18 G

Length 10 - 25 cm