

CHIBA NEEDLE

Can be considered as a Universal Needle used for Percutaneous Nephrostomy procedures. Cytological biopsies needle aspirations thoracentesis exploratory punctures to inject contrast dyes.



Key features:

- Enables soft penetration and avoids the risk of vessel laceration and reaches deep points without any excessive deviation.
- The special echogenic tip increases the visibility under the ultrasound guidance and guarantees the needles point localization
- Bevel tip needle.

Specifications

Size

Diameter 18-22-23 G

Length 15 - 20 - 22
cm