


## **COLLINS KNIFE - BIPOLAR**

Collin's Knife is used to make straight incisions during TURP Surgery for treatment of BPH. Product supplied Sterile and single use only.



### **Key features:**

- Faster firing, longer lasting
- Excellent performance
- High quality insulation.
- Compatible with all standard resectoscopes

### **Classification:**

Compatible with

- STORZ 1 Stem Instrument
- STORZ 2 Stem Instrument
- Olympus 1 Stem Instrument
- Wolf 1 Stem Instrument

## Other properties:

Available with / without bipolar cable.

## Specifications

### Size

**Diameter** 24 Fr