

COLLINS KNIFE - BIPOLAR

Collin's Knife is used to make straight incisions during TURP Surgery for treatment of BPH. Product supplied Sterile and single use only.



Key features:

- Faster firing, longer lasting
- Excellent performance
- High quality insulation.
- Compatible with all standard resectoscopes

Classification:

Compatible with

- STORZ 1 Stem Instrument
- STORZ 2 Stem Instrument
- Olympus 1 Stem Instrument
- Wolf 1 Stem Instrument

Other properties:

Available with / without bipolar cable.

Specifications

Size

Diameter 24

Fr