


## CUTTING LOOP - BIPOLAR

Bipolar Cutting Loop is used to remove prostate tissue in the form of small chips during TURP surgery for treatment of BPH. Product supplied Sterile and single use only.



### Key features:

- Faster firing, longer lasting.
- Excellent performance.
- High quality insulation.
- Compatible with all standard resectoscopes.

### Classification:

Compatible with

- STORZ 1 Stem Instrument
- STORZ 2 Stem Instrument
- Olympus 1 Stem Instrument

- Wolf 1 Stem Instrument

## **Other properties:**

Available with / without bipolar cable. Cutting Loops available in 0, 16, 20, 30, 45 and 90 Degrees.

## **Specifications**

Size

**Diameter** 24 Fr