


URETERAL STENTS OR DJ STENTS - SILICONISED

Used for internal drainage from the Uretropelvic junction to the bladder. Designed for long term indwelling.



Key features:

- Silicon - Polyurethane blended biocompatible material for patient safety and comfort
- Long term indwelling stents
- Marker bands for easy placements
- Highly Radiopaque for clear visualization
- Wide lumen ensures a superior drainage
- Stents are tapered at both ends, that helps in easy penetration to the exact location
- Can be kept indwelling for 365 days

Classification:

- Open End / Closed End
- Both End Multiloop / One End Multiloop
- Stent with Suture and / or Guide wire
- With / Without Hydrophilic coating

Specifications

Size

Diameter 4.5, 4.8, 5, 5.5, 6, 6.5, 7 and 8 Fr

Length 6 - 30 cm