

# DUAL LUMEN CATHETER

Used to insert guidewire and inject contrast simultaneously. It is used for cystometry and urethral pressure profilometry. Dual lumen catheterization and simultaneous bladder filling and pressure measurement. For single use only.



## Key features:

- Dual lumen provides an independent guide wire lumen and injection lumen
- The tapered tip enables the dilation of ureter
- Smooth surfaces for ease of introduction
- Eliminates the need for multiple catheterisation

## Specifications

Size

**Diameter** 8 -10 - 12

Fr

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email [info@acomedicaldevices.com](mailto:info@acomedicaldevices.com) [acomedicaldevices.com](http://acomedicaldevices.com)

**Length** 50 cm