


ESOPHAGEAL BALLOON DILATOR



Key features:

- Optimum radial force and uniform inflated balloon diameter.
- Radiopaque bands to identify balloon location fluoroscopically
- Atraumatic tip prevents damage to the endoscope
- Accepts 0.035" guidewire

Specifications

Size

Sheath Diameter 6.8 Fr

Balloon Inflated Diameter 6,8,9,10,11,12,13,15
mm

Length 1800 and 2400 mm

Balloon Length 55 mm

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acomedicaldevices.com acomedicaldevices.com

