

# ESOPHAGEAL DILATOR

Used to Dilate Esophageal and Colorectal strictures over previously placed wire guide.



## Key features:

- Soft and Flexible material for better precision and control
- Gradual Tapering of dilator for patient comfort
- Highly radiopaque distal and proximal markers
- Set supplied in a stylish and safe blow moulded carrying case
- Set includes dilators, stainless steel guidewire, cleaning brush 120 cm
- Available in set of 7 and 10 dilators
- Individual dilators provided in a peel open pouch
- Compatible with Ace Stainless Steel guidewire – 0.038 inch and 250 cm long

## **Classification:**

Set of 7 dilators. - Any 7 dilators from 5 mm to 20 mm

Set of 10 dilators - Any 10 dilators from 5 mm to 20 mm

## **Specifications**

### *Size*

**Diameter** 5 - 20 mm

**Length** 70 and 100  
cm