



ULTRA HD EUS- FNA NEEDLE

An Endoscopic Ultra Sound (EUS) - Fine Needle Aspiration (FNA) is used to sample targeted submucosal Gastrointestinal lesions through the accessory channel of an ultrasound endoscope. This device is supplied sterile and is intended for single use only. To be used by or on order of a physician only.



Key features:

- Enhanced Echogenic needle tip for optimal ultrasonic visualization.
- Unique needle design for higher and consistent tissue yield.
- Ergonomically contoured handle for better control and stability.
- Actuating lock for safer depth control of needle.
- Ultra smooth sheath with adjustable length attuned for a full range of endoscopes.

Classification:

EUS Needle Ordering Information

Needle Type	Needle			
	Size	19G	22G	25G
	Ø			

ACE ULTRA HD

FNA

AEUS-FNA-19

AEUS-FNA-22 AEUS-FNA-25

Specifications

Size

Sheath Diameter 5.2 Fr

Needle Size 19, 22 and 25 G