


GUIDE WIRES - STRIPED

Ace Striped Guidewires are used for percutaneous vessel entry to facilitate subsequent introduction(s) of an intravascular device. This device is Supplied sterile in a peel-open pouch. Intended for single use only..



Key features:

- Straight, Flexible tip - for atraumatic advancement through urinary tract
- Nitinol core with striped jacket - easy tracking of stents and urological instruments
- Low coefficient of friction for smooth advancement
- Better pushability for increased placement control
- Extremely smooth.
- Safe vessel access.

Classification:

with / without hydrophilic floppy tip

Specifications

E - 1/7 Neptune Society, Kalyani Nagar, Pune - 411006, Maharashtra, India. Tel +91 20 26681483, Fax +91 20 26681500

Email info@acemedicaldevices.com acemedicaldevices.com

Size

Diameter 0.018", 0.025", 0.028", 0.032", 0.035" &
0.038"

Length 150 cm