


## STONE BASKET - HELICAL

Helical Stone Baskets are designed for the manipulation and retrieval of stones from the Bladder, Ureter and Renal Pelvis. Intended for single use only.



### Key features:

- Four and five wire configurations are available
- Adjustable knob allows for controlled and assured basket expansion
- Handle permits complete dis-assembly
- Available in Stainless Steel and Nitinol
- Available in – PTFE, PEEK and Braided to suit every budget
- Choice of handles provided

### Classification:

- With Tip
- Tipless

## Specifications

### Size

**Diameter** 1.9 - 5 Fr

**Length** 70 - 130 cm